

# OWNER'S MANUAL

## ELECTRIC STERILMATIC STERILIZER

---

### DIGITAL VERSION

---

DOMESTIC MODEL:  STM-ED

EXPORT MODEL:  STM-EDX



- INSTALLATION
- OPERATION
- MAINTENANCE
- TROUBLE-SHOOTING
- PARTS & SERVICE
- WARRANTY



FORM NO.: **H-2117**

03/14 REV. A



Sterilizers.com - Your Sterilizer Expert  
Distribution - Maintenance - Guaranteed  
Alfa Medical - 10 Bond St Great Neck NY 11021  
1-800-762-1586 - info@Sterilizers.com

# **TABLE OF CONTENTS**

## **INTRODUCTION**

Product Description .....	i
Service Contacts .....	i

## **SECTION 1 INSTALLATION INSTRUCTIONS**

Service & Technical Information Contact .....	1
Operating Environmental Conditions .....	1
Intended Use - Sterilization Cycle .....	1
Installation .....	1
Electrical .....	1
Outside Venting .....	1
Water-Cooled Exhaust Condenser .....	1
Tray Supports .....	1
Baffle Installation .....	1
Operation Check .....	1
Cold Water Condenser .....	2
Service Connections .....	3
Electrical Requirements .....	3
Dimensions .....	3
Instructions for Installing Pan Supports & Baffles .....	4
Electrical Supply Connections .....	4
Wiring Diagram - 1 Phase .....	5
Wiring Diagram - 3 Phase .....	6
Typical Circuit Connection for STM-ED .....	7

## **SECTION 2 WATER CONDITIONS**

Water Conditions .....	9
------------------------	---

## **SECTION 3 OPERATING INSTRUCTIONS**

General Operating Instructions .....	10
Load Sterilizer .....	10
Close Door .....	10
Determine Correct Sterilization Times .....	10
Sterilization Guide .....	10
Minimum Sterilization Times Table .....	11

## **CONTROL PANEL**

Detailed Operating Instructions .....	12
Controller Overview .....	13
Set Up .....	13
Manual Programming .....	13
Setting Units .....	13
Setting Temperature .....	13
Setting Time .....	14
Setting Vent Mode .....	14
Stopping the Unit .....	14
Description of Cycles .....	14
Preset Keys .....	15
Programing Preset Keys .....	15
Using Preset Keys .....	15
Printing Data Logger .....	15
Print Key .....	16
Printing Previously Run Cycles .....	16
Printer Paper Cutting .....	16
Printer paper Changing .....	16
Low Water Reset .....	17
Steps to Reset the Unit .....	17
Recommended Steps but not Required .....	17

## **SECTION 4 MAINTENANCE**

Daily Cleaning Procedure .....	18
--------------------------------	----

## **SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY**

Sterilmatic Open Stand .....	19
Sterilmatic Open Stand with Condenser .....	19
Parts List for Condenser with Optional Stand .....	19
Door Adjustment .....	19
Door Assembly .....	20
The Door Gasket .....	20
Door Lift Spring .....	20
To Remove the Door Assembly .....	20
The Fulcrum & Drain Assembly .....	21
Roller Assembly .....	21
Cast-In Heating Elements .....	22
Low Water Cut-Off Manual Reset .....	22
Safety Valve .....	22
The Flue .....	22
Echhaust Solenoid Valve .....	22

## **SECTION 6 ILLUSTRATED PARTS LIST**

Vent Piping .....	23
Electrical Components .....	24
Sterilizer Assembly .....	25
Pan Supports & Baffle .....	26
Door Handle Assembly .....	27

## **SECTION 7 TROUBLE-SHOOTING**

Trouble Shooting Guide .....	29
------------------------------	----

## **SECTION 8 APPENDIX**

Error Codes .....	30
-------------------	----

## **SECTION 9 WARRANTY INFORMATION**

Sterilizer Warranty .....	31
Sterilizer Warranty Resistration .....	32

## **INTRODUCTION**

### **A. Product Description:**

The Market Forge Sterilmatic Sterilizer with Digital Controller (*model STM-ED*) is a compact, automatic, low cost steam pressure sterilizer (*autoclave*).

The sterilizing cylinder is a 3/16" (4.8mm) thick wall, welded aluminum. The exterior is made of polished stainless steel. Interior dimensions of 16" (406mm) diameter and 26" (660mm) long with a cubic content of 5,220 cubic inches (0.085 cubic meters) and has a door opening of 13-1/2" (343mm) wide and 11" (279mm) high. The sterilizing compartment has a pan capacity of:

- (3) 12" x 20" x 2 1/2" (305mm x 508mm x 64mm) or,
- (2) 12" x 20" x 4" (305mm x 508mm x 102mm) or,
- (1) 12" x 20" x 6" (305mm x 508mm x 152mm)

The sterilizer door is a self-sealing type that cannot be opened until the steam pressure is completely exhausted from the chamber. The door is 12 gauge stainless steel and removable for cleaning without tools. The door gasket is one-piece molded, also replaceable without tools or cement.

The sterilizing cycle is fully automatic with the time, temperature and venting controlled by the microprocessor based, digital controller.

The sterilizing temperature can be set anywhere in a range from 225°F (107°C) and 250°F (121°C). There is an on-board data logger/printer. The data logger records the time, temperature and pressure for each sterilization cycle. This data can be stored for future printing or printed out following each cycle.

### **B. Service Contacts:**

Should repairs be required, a network of authorized agencies is available to assist with prompt service. A current Directory of Authorized Service Agencies may be obtained by contacting:

Sterilizers.com - Alfa Medical - 10 Bond St Great Neck NY 11021  
1-800-762-1586 - info@Sterilizers.com

## **SECTION 1 INSTALLATION INSTRUCTIONS**

### **AUTOMATIC STERILMATIC STEAM PRESSURE STERILIZER MODELS: STM-ED, & STM-EDX**

#### **SERVICE & TECHNICAL INFORMATION CONTACT**

**NOTE:** *This unit should be serviced by qualified service personnel only.*

Your Sterilmatic Sterilizer has been developed to answer the need for a compact, automatic, low-cost steam pressure sterilizer. The following instructions cover installation. Should service be required, it is readily available by contacting our authorized service agency located nearest to you. The name of your local service company can be obtained by contacting the Service Department at Market Forge., Tel (617) 387-4100 or e-mail [custserv@mfii.com](mailto:custserv@mfii.com), or go to: <http://www.mfii.com/company/service> to find an authorized service agent in your area.

#### **OPERATING ENVIRONMENTAL CONDITIONS**

This unit is designed for commercial use and to be safe at least under the following conditions:

- For indoor use only.
- For use at altitudes up to 6500ft (2000m)
- For use at temperatures from 41°F (5°C) to 104°F (40°C).
- Maximum relative humidity 80% for temperatures up to 88°F (31°C) decreasing linearly to 50% relative humidity at 104°F (40°C).
- Main supply voltage fluctuations not to exceed ± 10% of nominal voltage.
- Transient overvoltages according to Installation Categories II (*in accordance with IEC 664*).
- Pollution Degree 2 (*in accordance with IEC 664*).

#### **INTENDED USE - STERILIZATION CYCLE**

This unit is intended to be operated intermittently. After a pre-heat cycle, the longest period of sterilization (heating) should be a maximum of 60 minutes. The digital timer allows up to 99 minutes, but it should be kept at no more than 60 minutes. After each use the unit should be opened for removal and reloading of product. The water level should be checked after each use and refilled when necessary.

#### **INSTALLATION**

Set sterilizer on counter, using the 6" (152mm) legs provided or assemble the optional stainless steel stand with under-shelf. If your Sterilmatic includes a water-cooled exhaust condenser, we recommend the use of the Sterilmatic stand, part number 95-6060. First, 1

level unit in place, then adjust rear legs to pitch the unit forward 1/4" (6mm) to insure positive drainage of the cylinder.

#### **ELECTRICAL**

Connect to proper electrical supply box and disconnect switch as shown on one of the following schematic diagrams - 208 or 240 volts, single or three phase. Connection is made from the rear of the unit, through the conduit to the terminal box located at the front of the unit. See installation specifications on page 5.

#### **OUTSIDE VENTING**

Connect 1/2" (13mm) nominal tubing exhaust to outside vent connection located on top of unit, within the control housing. **IMPORTANT:** Exhaust line must be vented to the outside to eliminate the exhausted steam and the accompanying noise from entering the room. Use 1/2" (13mm) copper tubing or suitable alternate. The overall height and length of the line should not rise more than 4 feet (1.2 meters) above the unit and exceed 15 feet (4.5 meters) with a minimum of bends. The line should slope downward after leaving the sterilizer in order to insure condensate drainage.

#### **WATER-COOLED EXHAUST CONDENSER**

If outside venting is not possible, an optional water-cooled condenser is available for connection to an open drain. If required order part no. 95-0436 kit.

#### **TRAY SUPPORTS**

Install side tray supports. Tray supports are attached by means of key-hole clearance slots which are slipped over studs located on the sides of the Sterilmatic chamber.

#### **BAFFLE INSTALLATION**

To insure maximum drying of packs, a baffle is supplied with your Sterilmatic. Place perforated splash baffle in bottom of the sterilizing chamber. Install small baffle with no perforation at the rear of the upper tray support channel.

**NOTE:** *The perforated baffle is not to be used as a shelf to place media or other items. It is intended to eliminate splashing.*

#### **OPERATION CHECK**

To check for proper operation of unit:

1. Close drain valve by turning handle clockwise.  
**WARNING:** DO NOT OPEN DRAIN VALVE WHILE UNIT IS OPERATING. PREMATURE OPENING MAY RESULT IN SCALDING OF OPERATOR.
2. Fill chamber with 4 to 5 quarts (3.8 to 4.7 liters) of ordinary tap water. DO NOT USE DISTILLED OR

## SECTION 1 INSTALLATION INSTRUCTIONS


DEIONIZED WATER.

3. Close chamber door.
4. Set exhaust selector to INSTRUMENTS AND PACKS (*fast exhaust*) or LIQUIDS (*slow exhaust*).
5. Set timer to 15 minutes. Cycle will go to completion automatically.

**NOTE:** Sterilizing cycle timer will not start until sterilizing set point temperature is obtained.

AMPS PER WIRE *STERILIZER					
Phase	3 Phase		1 Phase		208V (197-219) or 240V (220-240)
Volts	208	240	208	240	
Amps	26	30	45	52	

WIRING OF CONDENSING UNIT



**NOTES:**

1. Unit must be grounded and all wiring to comply with local codes.
2. Pipe to open drain. Do not make solid connection to sewer.
3. Condensing unit to be installed. Requires a front to back pitch, shown in the illustration below (see A).


### COLD WATER CONDENSER

ITEM	PART NO.	DESCRIPTION	QTY.
1	95-2119	Steam condensing unit	1
2	95-2219	Thermostat Box Assy.	1
3	95-0086	Exhaust line	1
5	15-7057	Copper tubing 3/8 OD	22.25"
6	10-1775	Rd. Hd. Mach. Screw, 1/4-20	2
7	10-2500	Lockwasher, 1/4	2
8	10-2308	Hex Nut, 1/4-20	2
9	10-1812	Rd. Hd. Mach. Screw, 10-32	2
10	10-2505	Lockwasher, 10	2
11	10-2340	Hex Nut, 10-32	2
12	95-4009	Front Template (7" Lg)	1
13	95-4010	Back Template (11" Lg)	1


### SERVICE CONNECTIONS REQUIRED

A	1/2" IPS Cold Water Connection
B	1" IPS Drain Connection (See Note 3)
C	115V Elec. Connection 7/8 Ø knockout (cond. unit)
D	Electrical Connection

## INSTALLATION INSTRUCTIONS - COLD WATER CONDENSING UNIT (KIT PART NO. 95-0436)



# SECTION 1 INSTALLATION INSTRUCTIONS




## REQUIRED CONNECTIONS:

A	Drain - 1/2" (13mm) FPT of 5/8" (16mm) OD copper (see note 1)
B	Steam Exhaust Connection - 3/8" (10mm) IPS (see note 2)
C	Electrical Connection - (*see table below)
D	Power Supply

**ELECTRICAL REQUIREMENTS:** Unit will be rated at 12kW at 236V and will be equipped for operation at:

- 208-240V, 3 Wire, 1 Phase, 60 Hz
- 208-240V, 4 Wire, 3 Phase, 60 Hz
- 220V, 1 Phase, 50 Hz
- 220V, 3 Phase, 50 Hz



## AMP/PHASE

DOMESTIC MODEL	kW	Hz	Single Phase		Three Phase	
			208V	240V	208V	240V
STM-ED	9.3	60	45A	-	26A	-
	12.4	60	-	52A	-	30A

EXPORT MODEL	kW	Hz	Single Phase		Three Phase	
			220V	240V	220/380V	240/415V
STM-EDX	10.4	50	48A	-	16A	-
	12.4	50	-	52A	-	18A

- **UNIT MUST BE GROUNDED**
- **MAIN SUPPLY VOLTAGE FLUCTUATIONS ARE NOT TO EXCEED ± 10% NOMINAL SUPPLY VOLTAGE**


## NOTES:

1. An air break must be provided if a unit drain line is run.
2. Vent exhaust to atmosphere. B1 is actual connection, but must exit casing at B.

**IMPORTANT:** Exhaust line must be vented to the outside to eliminate the exhaust steam and the accompanying noise from entering the room. Use 1/2" (13mm) copper tubing or suitable alternate. The overall height and length of the line should not rise more than 4' (1.2 meters) above the unit and exceed 15' (4.5 meters) with a minimum of bends. The line should slope downward after leaving the sterilizer in order to ensure condensate drainage.

**IMPORTANT:** Failure to comply with this outline will affect the sterilization process.



When an exhaust condenser is supplied; the following services must be provided: 1/2" (13mm) IPS cold water: 1" (25mm) IPS waste: 115V electrical line.





## SECTION 1 INSTALLATION INSTRUCTIONS

### INSTRUCTIONS FOR INSTALLING PAN SUPPORTS AND BAFFLES

1. Locate the mounting studs on the inside of the chamber. There are two rack mounting studs on each side (see Fig. 1).
2. Taking one pan support and positioning rack so that the pan stop is facing the rear of the unit and the wires are facing toward the center of the unit. The pan stop is a piece of sheet metal welded to the rack with a 65° bend.
3. Begin to hang the pan support by placing the rear key-way slot onto the rear mounting stud and slide the rack until the slot sits on the mounting stud. When this is done correctly the front mounting stud will be in position to place the front key-way slot. Slide the rack down into its correct position (see Fig 2).



4. After installing one pan support rack correctly, you can install the upper baffle. Position the baffle so that the 45° bend is facing up towards the front of the unit (see Fig. 3). Slide the mounting tab onto the flat bend on the pan stop bracket. The baffle should now stay in place by itself, but in a tilted state (see Fig. 4).
5. Position the second pan support rack into the cavity and slide the other mounting tab onto the rack flat bend while the pan support rack is not on the mounting studs. Hang the pan support by placing the rear key-way slot onto the rear mounting stud and slide the rack until the slot sits on the mounting stud. When this is done correctly the front mounting stud will be in position to place the front key-way slot. Slide the rack down into its correct position.



6. Place the Perforated Water Baffle so that it sits on the bottom of the inside of the sterilizer chamber (see Fig. 5).
7. **CAUTION:** Do not cover the holes in the Perforated Water Baffle by using it as a shelf. This will result in a disrupted flow of steam.

**THE ELECTRIC SUPPLY CONNECTIONS FOR STM-ED:** Connect to proper electrical supply as indicated on nameplate on top of unit. The power supply cord is brought in from the rear of the unit, through the conduit and the connection is made at the terminal box located at the front of the unit.


**THE ELECTRIC SUPPLY CONNECTIONS FOR STM-EDX:** Connect to proper electrical supply as indicated on nameplate on top of unit. Connection is made from the rear of the unit, through the conduit to the terminal box located at the front of the unit. All control circuits are 220 volts.

In order to accomplish this, a current-carrying grounded neutral must be provided.

Thus, a three phase system must be 4-wires. Most electrical codes require, and we recommend, that a separate switch be located within sight of the sterilizer.

# SECTION 1 INSTALLATION INSTRUCTIONS

Volts	Phase	AMPs Per Line Wire			kW
		L1	L2	L3	
208	1	45	45	--	9.3
208	3	26	26	26	9.3
220	1	48	48	--	10.4
220	3	28	28	28	10.4
240	1	52	52	--	12.4
240	3	30	30	30	12.4




DIGITAL STERILMATIC MODEL STM-ED - 1 PHASE WIRING DIAGRAM  
DRAWING No. 95-6311 Revision A - Date 10/25/13

DRAWING No. 95-6311 Revision A



# SECTION 1 INSTALLATION INSTRUCTIONS

Volts	Phase	AMPs Per Line Wire			kW
		L1	L2	L3	
208	1	45	45	--	9.3
208	3	26	26	26	9.3
220	1	48	48	--	10.4
220	3	28	28	28	10.4
240	1	52	52	--	12.4
240	3	30	30	30	12.4




DIGITAL STERILMATIC MODEL STM-ED - 3 PHASE WIRING DIAGRAM  
 DRAWING No. 95-6312 Revision A - Date 10/25/13

DRAWING No. 95-6312 Revision A

# SECTION 1 INSTALLATION INSTRUCTIONS

## TYPICAL CIRCUIT CONNECTION FOR STM-ED STERILIZER



## **SECTION 1 INSTALLATION INSTRUCTIONS**

This page intentionally left blank.

## **SECTION 2 WATER CONDITIONS**

Market Forge from time to time is asked the question about using distilled or deionized water for use with our Sterilizer models STM-ED and STM-EDX. We are always asked why these water choices are not allowed for use with our units and what would be recommended. To address this situation, we have compiled the following as a means of satisfying these questions:

1. We have found that the use of distilled or deionized water will aggressively attack the pure coat of Aluminum Alclad, which protects the bottom surface from oxidizing and then eventually pitting (*reference: Operating and Maintenance Instructions*).
2. In addition pitting can also be caused by several other external environmental factors. Few examples are as follows. These conditions have been highlighted in our documentation.
  - Grains of hardness in the water supply should be as follows (.25 to 2).
  - A pH imbalance in the water supply can greatly affect the life to the aluminum cylinder. The pH range that would be recommended is between 7.0-8.5.
  - The lack of a positive electrical ground can cause an electrolytic reaction that will accelerate pitting.
  - Another contribution to accelerate pitting is the type of cleaning solutions used or the abrasive scrubbing pads. If a low pH is present with the detergents being used or an abrasive pad, the protective Alclad coating will be removed during the cleaning process.
  - Spillage of media being sterilized can also contribute to the accelerated pitting if it is corrosive.
  - CHLORINE LEVEL  $\leq$  1 PPM.

### **IMPORTANT NOTE:**

***Market Forge will not be responsible for damage resulting from the use of hard or corrosive water, from failure to drain the unit daily, or from inadequate cleaning procedures.***

## SECTION 3 OPERATING INSTRUCTIONS

### GENERAL OPERATING INSTRUCTIONS:

- 1. IMPORTANT:** Make sure the drain valve is closed. Fill bottom of the sterilizer chamber with approximately six quarts of water or just below ledge at bottom of door opening. (*If water supply is known to be hard or corrosive, a source of treated water should be used.*) DO NOT USE DISTILLED OR DEIONIZED WATER. (*See section 2*)
- 2. LOAD STERILIZER:** Use proper sterilizer loading procedures when placing materials in sterilizer chamber. All solid containers or instruments must be placed so that water or air will not be trapped in them.
- 3. CLOSE DOOR:** Grasp handle, and holding it in vertical position, pull door down until bottom of door rests in the bottom of door opening. Then rotate handle forward, engaging the lower curved portion under the horizontal rod in the casting at the bottom of the door opening. Push handle all the way down and back until door is locked securely in position.
- 4. DETERMINE CORRECT STERILIZATION TIMES:** (*Referring to page 11 for minimum sterilization times table.*) **NOTE:** *In no case should the timer be set to less than 15 minutes. Sterilization will not be accomplished in less than 15 minutes exposure time.*
- When the sterilizer chamber reaches the selected temperature, the timed Heating/Sterilization cycle will begin. When the Heating/Sterilization cycle is completed, the electric supply to the heating elements will be opened (shut off) automatically. When the chamber pressure reaches 0 (zero) the door may be opened.  
**NOTE:** *Before opening the chamber door be sure to have the Control Panel Flip Cover in the 'DOWN' position. This protects the LCD screen from coming into contact with too much steam.*  
  
At this point you may release the handle and let go to avoid possible contact with the remaining escaping steam. When opening the door allow a few seconds for steam to escape from the chamber before opening completely.  
  
**NOTE:** *For more detailed Operating Instructions please refer to Section 3 titled 'DETAILED OPERATING INSTRUCTIONS'*
- To assist in drying racks, release door handle after pressure has been attained at start of cycle. Pressure in chamber will keep door closed. The use of a wire basket will provide better drying for dressings. At end of sterilizing cycle, release door handle and open slightly. Do not lift door to open

position. This will allow steam and moisture to escape. Allow door to remain in this position for 15 to 20 minutes before removing load. Small packs can be dried successfully with this procedure. We do not recommend the sterilization of large packs, such as linens. Be sure condensate baffles are in position in the chamber.

7. Remove load and check water level for next operation.

### STERILIZATION GUIDE:

- **PACKS (*Linens, gloves, etc.*):** Use wire basket to facilitate drying. Be sure condensate baffles are in place. Place packs on edge and arrange load in chamber, so that only minimal resistant to passage of steam through the load will exist.

**NOTE:** *Place gloves in upper two-thirds of chamber.*

- **JARS, CANISTERS (*etc.*):** Place containers on side to allow for displacement of air and complete contact of steam to surfaces. Drying is also facilitated.
- **PETRI DISHES, PIPETTES, DESICCATORS (*etc.*):** Should be inverted.
- **UTENSILS, TREATMENT TRAYS:** Placed on edges to facilitate drying.
- **INSTRUMENT SETS:** Place instruments set in trays having mesh or perforated bottoms. Place trays flat on shelves.
- **COMBINING FABRICS & HARD GOODS:** Place hard goods on lowest shelves.
- **PLASTIC UTENSILS:** DO NOT stack or nest plastic items.
- **LIQUIDS:** Sterilize liquids separately from other supplies or materials. Set vent to slow.
- **SMALL ITEMS:** Sterilize small items in baskets, or trays.

**NOTE:** *IF THE EQUIPMENT IS USED IN A MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.*

## SECTION 3 OPERATING INSTRUCTIONS

### MINIMUM STERILIZATION TIMES


TIME ( <i>Minutes</i> )	ARTICLES
15	<ul style="list-style-type: none"><li>• Glassware, empty, inverted.</li><li>• Instruments, metal in covered or open tray, padded or unpadded.</li><li>• Needles, unwrapped.</li><li>• Pipettes, blood diluting, serological, volumetric, etc</li><li>• Tubing glass (<i>6mm</i>), (<i>10mm</i>) inverted</li></ul>
20	<ul style="list-style-type: none"><li>• Flasked solutions 75-250 ml.</li><li>• Instruments, metal combined with other materials in covered and/or padded tray.</li><li>• Instruments wrapped in double thickness muslin.</li><li>• Rubber gloves, catheters, drains, tubing, etc. Unwrapped or wrapped in muslin or paper.</li></ul>
30	<ul style="list-style-type: none"><li>• Brushes in dispensers, in cans of individually wrapped.</li><li>• Dressings, wrapped in paper or muslin, small packs only.</li><li>• Flasked solutions 500-1000 ml.</li><li>• Syringes, unassembled, individually packaged in muslin or paper.</li><li>• Needles, luer, individually packaged in glass tubes or paper.</li></ul>
45	<ul style="list-style-type: none"><li>• Flasked solutions 1500-2000 ml.</li></ul>













## SECTION 3 OPERATING INSTRUCTIONS

### DETAILED OPERATING INSTRUCTIONS

Before operating this unit be sure to have read the Owner's manual for proper setup, service connections and installation. In addition, the Owner's manual will cover sterilization recommendations, daily cleaning procedures and parts lists.

#### A. Control Panel Identification



ITEM	DESCRIPTION	BUTTON
1	Control Panel Power - ON ( 1 ), OFF ( 0 )	
2	Reset Button	
3	Low Water Indicator Light	
4	Down Key	
5	Printer Key	
6	Preset Keys	
7	Up Key	
8	Start / Cancel Key	
9	Digital LCD Display	
10	Temperature Key	
11	Time Key	
12	Vent Key	
13	Print out	
14	Printer Door Button	

## **SECTION 3 OPERATING INSTRUCTIONS**

### **B. Controller Overview**

The Digital Controller is made up of the Operational Keys, Digital LCD Display and the Data Logger/Printer. The operator can set the controller to display in either degrees Fahrenheit ( $^{\circ}F$ ) with pressure in PSI or in degrees Celsius ( $^{\circ}C$ ) with pressure in kPa.

During operation the Digital LCD Display will show the actual temperature within the sterilizer chamber, the time remaining for the sterilizer cycle and the vessel pressure in digital form and bar form. It will also display the current cycle state, Heating up, Timing (*sterilizing*) or Venting.

The three operating parameters that can be set by the operator are sterilizing **Temperature**, sterilizing **Time** and **Venting** mode.

The sterilizing temperature can be set in a range from 225 $^{\circ}F$  (107 $^{\circ}C$ ) to 250 $^{\circ}F$  (121 $^{\circ}C$ ).

The sterilizing Time can be set in Minutes and Seconds (*if desired*).

The Venting mode can be set to **FAST** (*approximately 3 minutes*) or **SLOW** (*approximately 11 minutes*).

To initiate a sterilizing cycle you need to set three parameters, **Temp**, **Time** and **Venting**. Or you can set commonly used values for these parameters into three 'Preset' keys. Each Preset key can store a value for all three parameters. The parameter values stored in the Preset keys can easily be selected before running a sterilization cycle. This avoids you having to set all three parameters every time you run the unit.

Each sterilization cycle records the temperature, time and pressure at one minute intervals during the heat-up cycle, the sterilization cycle and the venting cycle. These cycles are described in more detail below. The recorded data can be printed out after each complete run or printed out later, on an as-needed basis.

### **C. Setup**

Make sure unit is connected to its electrical source. Make sure unit has recommended amount of water in its chamber. Make sure unit drain is closed. If a ventilation hood is required, make sure unit is placed accordingly under the hood.

Be sure that the unit's door is closed and locked before operating!

Once the unit has power, push the black Control Panel Power switch to On ( I ). This will bring power to the control panel and illuminate the Digital LCD Display screen.

You are now ready to use the unit.

### **D. Manual Programming**

The first step is to decide what form of temperature and pressure measurement you want to use, Fahrenheit and PSI or Celsius and kPa. This can be changed at any point during operation.

- **Setting units** ( $^{\circ}F$  & PSI) or ( $^{\circ}C$  & kPa);
  - Press and hold the **UP** key and the **DOWN** key while simultaneously pressing and releasing the **TEMP** key will toggle between ( $^{\circ}F$ /PSI) and ( $^{\circ}C$ /kPa).
- **Setting Temperature**
  - While the unit is powered up but not in use, considered the 'Idle State', the temperature being displayed is the current temperature within the sterilizer chamber.
  - Pressing and releasing the **TEMP** key will display the current set-point temperature for three seconds.
  - To set the target sterilization temperature set-point;



## SECTION 3 OPERATING INSTRUCTIONS

- Press and hold the **TEMP** key until the display goes blank and shows only the temperature. Now press the appropriate **UP** or **DOWN** key to get to desired set-point.
- When the full display returns the target set-point has been set.
- **Setting Time**
  - While in the 'Idle State' the Time being displayed is the current time set-point.
  - To set the target sterilization time set-point;
    - Press and hold the **TIME** key until the display goes blank except for the minutes. Now press the appropriate **UP** or **DOWN** key to reach the desired sterilization set-time minutes.
    - After three seconds the minutes will blank and the seconds set-point will be displayed. If you need to set the desired seconds, again, use the appropriate **UP** or **DOWN** keys.
    - If seconds are not required than either don't touch any key and after three seconds the full display will come back and the **TIME** is now set or immediately press the **TIME** key again and that will also set the time.
- **Setting VENT mode**
  - While in the 'Idle State' the **VENT** mode will be displayed as either **FAST** or **SLOW**, whichever is active.
  - To set the desired **VENT** mode (*SLOW* or *FAST*);
    - At any point during any cycle the **VENT** mode can be toggled simply by pressing and releasing the **VENT** key.
- Starting the unit
  - Once you have set your desired **TEMP**, **TIME** and **VENT** type you can now proceed in starting the unit.
  - To start unit;
    - Press the **START/CANCEL** key
- **Stopping the unit**
  - Once the unit is started it will run through all 3 cycles automatically (*cycles are described below*) and finish in the **DONE** state (*also described below*).
  - If at any time throughout the cycle you have to stop or cancel the cycle just press and hold (*for 3 seconds*) the **START/CANCEL** key.
  - After stopping the unit in this manner the unit will still have to go through the **VENT** cycle to release the pressure within the chamber before you can open the unit door.
- **Description of cycles;**
  - **HEATING** cycle – Once the **START/CANCEL** key is pressed the unit will flash the **HEATING** icon while the sterilizer is heating up to the sterilization temperature set-point. While it is heating the set-point TIME is displayed and will not change.
  - **TIMING** (*sterilizing*) cycle – When the unit reaches its set-point temperature the unit will enter the sterilization or **TIMING** cycle. At this point the unit will display and count down the sterilization time.
  - **VENTING** cycle – When the timing (*sterilization*) cycle time reaches zero the set **VENTING** cycle begins. The **FAST VENT** is programmed to vent for 3 minutes through the solenoid valve. The **SLOW VENT** is programmed to vent for 11 minutes through a bleeder orifice.
  - **DONE** state – At the completion of the **VENTING** cycle the unit will display **DONE** and a beeper will sound. The beeper sequence will be **ON** for 1 second and **OFF** for 9 seconds. This sequence will repeat until any key on the control panel is pressed or 3 minutes have passed.

Once in the **DONE** state it is safe to open the chamber door. Before opening the chamber door be sure to have the Control Panel Flip Cover in the 'Down' position. This protects the LCD screen from coming into contact with too much escaping steam from the chamber. When the door is left in the closed position you may notice a rise in the temperature and pressure on the LCD screen. You should open the chamber

## SECTION 3 OPERATING INSTRUCTIONS

door to allow the chamber to cool.

**NOTE:** *There is never any harm in releasing the chamber door latch. If there is pressure inside the chamber the door cannot open due to its design which does not allow the door to open under pressure.*

### E. Preset Keys

There are three **PRESET** keys numbered **1, 2 & 3**. In each key you can save set-point values for the Temperature, Time and Vent mode. For example, if you commonly need to sterilize a media that requires specific set-point values, you can save these specific values into a **PRESET** key. This way you don't need to keep setting the Temp, Time and Venting info before each running of a sterilization cycle.

#### • Programming PRESET Keys

- To save values into a PRESET key;
  - Press and hold the desired **PRESET** key (1, 2 or 3). The display shows the values currently stored in that preset key. When the “**PRESET**” icon and the preset “number” icon in the display start to flash, the values are now ready to re-program. At this point you can change any or all of the parameters, Temp, Time, Vent as described previously.
  - Every three seconds without an **UP** or **DOWN** key press the display will proceed from **TEMP**, **MINUTES**, then to **SECONDS** waiting for new preset set-points then exits the program preset state.
  - During this programming state you can advance from the Temp menu to the Time menu simply by pressing the **TIME** key.
  - Likewise, while in the Time menu you can advance from Time-minutes to Time-seconds simply by pressing the **TIME** key.
  - In addition, while in the Time menu you can go back to the Temp menu by pressing the **TEMP** key.
  - At any point during this programming state you can press the appropriate **PRESET** number key and that will store the new values and exit the programming state.

#### • Using PRESET Keys

- While the unit is in the 'Idle State' simply press the appropriate **PRESET** key (1, 2 or 3) to use the values previously stored into that **PRESET** key. At this point you just need to press the **START/CANCEL** key.

### F. Printing/Data Logger

This unit is set up to record the Time, Temperature and Pressure during all phases of the full sterilization cycle. All three of these parameters are recorded at one-minute intervals. The data logging clearly labels the printed output with headers separating the three cycles, **HEATING**, **STERILIZING** and **VENTING**.

Each run is assigned a unique ID number that is printed at the top of the printout. Each unique ID number is made up of two numbers separated by a dash (-).

(Example; #237-00053)

**PRINT OUT EXAMPLE** ⇨

#B41-00003		
HEATING		
TIME	F	PSI
0:00	72	0.0
1:00	72	0.0
2:00	168	0.0
3:00	201	0.0
4:00	214	0.7
5:00	246	12.9
5:26	250	15.1
STERILIZING		
TIME	F	PSI
0:00	250	15.1
1:00	251	15.6
2:00	251	15.6
3:00	251	15.6
4:00	251	15.7
5:00	251	15.8
6:00	251	15.8
7:00	251	15.6
8:00	251	15.6
9:00	251	15.6
10:00	251	15.6
11:00	251	15.6
12:00	251	15.6
13:00	251	15.6
14:00	251	15.6
15:00	251	15.6
16:00	251	15.6
17:00	251	15.6
18:00	251	15.6
19:00	251	15.6
20:00	251	15.6
21:00	251	15.6
22:00	251	15.6
23:00	251	15.6
24:00	251	15.6
25:00	251	15.6
26:00	251	15.6
27:00	251	15.6
28:00	251	15.6
29:00	251	15.6
30:00	251	15.6
VENTING		
TIME	F	PSI
0:00	251	15.6
1:00	212	0.0
2:00	212	0.0
3:00	212	0.1
TOTAL		
38:26		
DATE: _____		
OPERATOR: _____		
COMMENTS: _____		

The first 3 digit number is a number that never changes. It is a unique number that identifies the sterilizer unit itself. The last five digits are incremented by one for every time the unit is run.

## SECTION 3 OPERATING INSTRUCTIONS

Also included at the end of each printout is the following;

Example: **DATE:** \_\_\_\_\_  
**OPERATOR:** \_\_\_\_\_  
**COMMENTS:** \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

The DATE, OPERATOR and COMMENTS fields can be hand written in if required.

- **PRINT Key**

- In the LCD Display, if the word PRINT is displayed then the recorded data will be automatically printed at the end of the full sterilization cycle. Likewise, if the word PRINT is not displayed then the recorded data will not be printed automatically after the full sterilization cycle completes.
- To toggle between PRINT enable and PRINT disable just press and release the PRINT key. This can be performed at any time during the Heating, Timing (*Sterilization*) or Venting cycles.

- **Printing previously run cycles**

- The data logger will store approximately 20 full sterilization cycles. This number may vary slightly depending on how long each cycle runs. For example, the data logger can store more cycles set at 30 minutes versus cycles all set at 60 minutes.
- As stated before each full cycle is assigned a unique ID number. This number can be used to scroll back through the data logger and print out the exact cycle you want to print out. There is also an option to print out ALL the previous cycles.
- To print out selected cycles;
  - Press and hold the PRINT key. The display will show the ID of the last run cycle. To scroll back through the cycles just press the DOWN key. Likewise, during the scrolling, you can hit the UP key to scroll forward through the cycles. At your selected cycle ID number just press and release the PRINT key.
- To print out all the cycles;
  - Press and hold the PRINT key. When the display shows the last run ID number keep hitting the DOWN key until the word 'ALL' is displayed. Now press and release the PRINT key.
  -

- **Printer Paper cutting**

- The paper will be cut automatically at the end of each single print or printing ALL.

### **G. Printer Paper Changing**

Follow these steps to re-load the printing paper;

1. Push the printer door "open" button. Two doors will open, the inner door and outer door.
2. Remove paper roll spool from printer. Remove black plastic spindle from center of spent paper roll. **DO NOT DISCARD BLACK PLASTIC SPINDLE!**
3. Insert black plastic spindle into new paper roll. Insert new paper roll into printer. (**Note: paper exits from top of roll**)
4. Hold paper tip with one hand and close the inner door by pressing on the yellow strip.
5. Close black outer door making sure that paper protrudes through the slot.



## **SECTION 3 OPERATING INSTRUCTIONS**

### **H. Low Water Reset**

If the water inside the chamber is allowed to run dry it will trigger the Low Water Cut-off. At this point the unit will;

- Shut down all three heating elements
- Light the red Low Water indicator light.
- Start the Venting cycle
- The Data Logger will record the error code

*(Note: The sterilizer's LCD display may show the temperature and pressure rising slightly during this venting period but that is normal.)*

- Once the unit completes the VENTING cycle the screen will flash 'DONE' and display 'SERVICE' and the Beeper will sound.

### **Steps to Reset the unit;**

1. The LCD display now shows 'ERR 06' and beeper changes to a constant tone.

*(Note: An explanation of Error Codes is shown below in the Appendix)*

### **Recommended steps but not required:**

- Shut off 'ON/OFF' power switch to controller.
  - Wait until unit cools down. Opening the unit door will help the unit to cool down quicker.
  - When power is restored the buzzer tone will continue and 'ERR 06' will be displayed again.
2. Open the door and add water to the chamber.
  3. Once the chamber has cooled enough you can press the Low Water Cutoff RESET button. Listen for the 'click' sound.
  4. If constant tone continues press any menu key to turn off.
  5. The unit is now back in the 'Idle State' and ready to be run again.

## **SECTION 4 MAINTENANCE**

### **DAILY CLEANING PROCEDURE (AT THE END OF EACH DAY):**

1. Remove bottom splash baffle.

***NOTE: IMPORTANT! STERILIZING CHAMBER MUST BE CLEANED AND DRAINED DAILY USING THE FOLLOWING PROCEDURE. WASH WETTED PORTION OF THE CYLINDER THOROUGHLY BY ADDING A MILD DETERGENT TO WATER IN CYLINDER.***

2. If a soft cloth or brush is used with the detergent and does not completely remove the surface film, a nylon soap pad should be used. After washing thoroughly rinse with clean water. Dry cylinder\* and leave door open overnight.

\* The Sterilmatic cylinder is constructed of corrosion resistant Alclad aluminum alloy. The protective properties of this material afforded to the interior portion of the cylinder which is exposed to water may be destroyed by allowing a film to form. Such a film can be caused by salts or other contaminants in the water. Corrosion may also occur if water is not drained daily.

## SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY



**STERILMATIC OPEN STAND:**  
Market Forge Sterilmatic Stand can be supplemented with an Optional Stand for utility use where maximum compactness is desired.

The sturdy, stainless steel unit is equipped with adjustable leg extensions which allow the unit to be installed and leveled over existing contours in the floor.

The open design lends itself to maximum sanitary conditions because of the ease with which periodic cleaning can be done.

Though simple in design and appearance, the sterilmatic stand is the ideal arrangement for mounting in that it allows secondary air to circulate.

### **STERILMATIC OPEN STAND WITH CONDENSER:**

Market Forge can provide the open stand with an optional steam condenser system for use where steam exhaustion into the room is undesirable.

The condenser is automatically controlled by the thermostat. The normal factory thermostat setting is 130°F (54°C). The open under-shelf of the stand gives added utility providing a handy tabouret for utensils and access for drainage of water from the sterilizing chamber.


### **PARTS LIST FOR CONDENSER WITH OPTIONAL STAND:**

PART NO.	50 Hz	DESCRIPTION
10-4653	10-4653	Thermostat
10-4035	10-7074	3/8" Solenoid
10-5731	10-5731	1/2" Water Stop Valve
95-2106	95-2106	Water Injection Assy.
95-1680	95-1680	Shelf

### **STERILMATIC DOOR ASSEMBLY:**


The Door of the Sterilmatic has been engineered to establish a positive method of sealing the steam pressure within the sterilizing cylinder. As steam pressure builds up within the cylinder, the door seal will tend to become more positive.

However, the door should be adjusted to make a good initial seal between the door gasket and the door opening without the added assistance of internal cylinder steam pressure with the simple action of securing the door handle down in a locked position, the door gasket should be sufficiently compressed against the door opening, all the way around to prevent any steam leakage from occurring.




### **DOOR ADJUSTMENT:**

The Door Adjustment is Located in the Fulcrum Casting at the base of the door opening. This adjustment employs the use of a screw and locknut in order to adjust the Sterilmatic Door to a tighter closed position (to prevent steam from leaking by the door gasket as pressure builds up), it is necessary to loosen the locknut and back off the screw at least one-quarter of a turn and re-tighten the locknut (see Door figure 1 shown above & Door Figure 2 shown below).



FOR DOOR SEAL ADJUSTMENT  
 • TO TIGHTEN DOOR •  
 LOOSEN 'A' AND TURN 'B'  
 1/4 TURN COUNTERCLOCKWISE.  
 • TO LOOSEN DOOR •  
 TURN 'B' 1/4 TURN CLOCKWISE.  
 BE SURE TO TIGHTEN LOCKNUT  
 'A' AFTER ADJUSTMENT.

## SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY



### DOOR ASSEMBLY

ITEM	PART NO.	DESCRIPTION
1	10-6765	Pivot Spring Bearing
2	91-2718	Right and Left Door Spring ( <i>sold as a pair</i> )
3	10-1776	10-32 Machine Screw 1/2" Long
4	10-2666	Door Gasket
5	95-3204	Door & Door Spring Assembly
-	95-0124	Items 1 through 5

#### THE DOOR GASKET:

Keep the gasket clean. With normal closing and locking of the door assembly, a steam-tight seal should be made between the door gasket and the door opening. This seal cannot be maintained if particles of foreign matter are allowed to accumulate upon either of the contacting surfaces.

If there is leakage by the door gasket before a steam build-up within the steam chamber and leakage does not stop when the sterilizer reaches sterilizing temperature and pressure than regard the door assembly as improperly adjusted. A re-adjustment must then be made of the seal adjustment door screw.

To change the door gasket, remove the entire door assembly as a unit. Discard the old gasket, replace it with a new one (no cement is required), and reinstall the door assembly. Make an operational check for leakage and adjust the door, if necessary.

#### DOOR LIFT SPRING:

Market Forge supplies door lift springs in sets only. This policy has been found to be in the best interest of the customer. Through continuous use, some of the original qualities of the springs are lost and it becomes advantageous to make replacements to both the left and right door lift springs in the event that one becomes damaged or broken.

Replacement door lift springs are marked with tabs at the factory prior to shipment to identify a right from a left spring. These springs must be installed with the right door lift spring on the right of the door and the left door lift spring on the left of the door as viewed from the front of the sterilizer.

#### TO REMOVE THE DOOR ASSEMBLY:


The Door Assembly can be removed from the inner sterilizing chamber as a unit without the use of any special tools or equipment. However, a systematic approach to this is warranted as the clearances through the portal are close, and much confusion can result if not removed in the sequence described below:

1. First, lift off and remove the two pan supports to expose the door linkage on either side of the inner sterilizing chamber .
2. Raise the door to a fully opened position, and disengage the door spring from each of the door spring studs. Accomplish this by counteracting the force of the door lift spring with one hand while working the end of the door spring off the spring stud with the free hand. Do this on both sides of the door assembly.
3. When the end of the door springs have been completely freed from their respective door spring studs, the door springs on either side of the door assembly can easily be slipped off their studs.
4. Rotate the entire door assembly out through the door opening, passing the door handle through the opening first, and then one end of the door spring as shown in the illustration. The remainder of the door assembly will then pass through the door opening quite easily.
5. To replace the door assembly, reverse the step-by-step procedure described above.

## SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY

### **THE FULCRUM & DRAIN ASSEMBLY:**

The fulcrum and drain assembly is located at the lower front of the sterilizing chamber and furnishes a sturdy anchorage for the door locking system of the door handle. Also provided in this assembly is a means for adjustment of the door seal. The drain port and drain valve provide a means of discharging accumulations of water from within the sterilizing chamber.



### **ROLLER ASSEMBLY (Items 8 & 9):**

The Roller Assembly must be kept free-rolling at all times. Should this assembly be allowed to become frozen due to lack of lubrication, undue strain will be put on the door handle and the fulcrum casting while the door is being locked. Use only a dry lubricant such as graphite; as oil or grease will tend to attract dirt to this area.

### **FULCRUM & DRAIN ASSEMBLY**

ITEM	PART NO.	DESCRIPTION
1	10-3116	1/4" - 20 X 5/8 helicoil
2	10-1999	10-32 Machine screw, 1 5/8" long
3	10-2358	1/4" - 20 fulcrum nut
4	10-2087	1/4" - 20 allen set screw
5	10-3111	1/4" - 20x 3/8 helicoil
6	10-2513	1/4" Shakeproof washer
7	10-1763	1/4" - 20 Machine screw 3/4" long
8	95-0120	Bearing spacer
9	95-0198	Bronze Bearing
10	10-3111	1/4" - 20 x 3/8 helicoil
11	10-2513	1/4" Shakeproof washer
12	10-1790	1/4" - 20 Cap screw 7/8" long
13	10-4485	Drain valve knob
14	10-2514	#10 Shakeproof lockwasher
15	10-2318	10-32 acorn nut
16	95-2643	Adapter - steinball valve
17	10-1950	6-32 Round head screw 1 5/8" long
18	95-2616	Front outer case lower
19	95-0116	Fulcrum and drain casting
20	10-1049	Nipple 1/2" IPS 2 1/4" long stainless steel
21	10-1041	Ball valve stein
-	95-0115	Fulcrum and drain assembly, Items 1 through 12, 14, 15, and 19



## **SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY**

### **CAST-IN HEATING ELEMENTS:**

Located under the sterilizing cylinder is a bank of (3) U-shaped heating elements. These elements are welded in place in a protective aluminum shield. The elements cannot be removed, and in the unlikely event that one or all fail, the complete cylinder must be replaced.

### **THE LOW WATER CUT-OFF (MANUAL RESET):**

Fastened to a special mounting brace behind the front panel, the Low Water Cut-Off acts to shut off the complete unit, should the water run dry. The Low Water Cut Off is factory set, to shut the unit off when the cylinder temperature rises between 380° and 440° Fahrenheit.

When the Sterilmatic is turned on without water or the water has been evaporated away, the temperature of the aluminum sterilizing cylinder will rise and by heat induction effect the Low Water Cut-Off. Its inner electrical contacts will be forced open from heat expansion, thus cutting off the flow of electric current to the heating elements. With the replacement of water into the cylinder the cylinder temperature will drop and the contacts of the Low Water Cut-Off can be again closed. The unit will only restart after the manual button has been re-set.


### **THE SAFETY VALVE:**

The Safety Valve is factory set to automatically open and exhaust excess steam from within the sterilizing cylinder, thereby assuring that operating pressures remain within safe limits. The lever action of the safety valve must be free to operate unrestricted at all times. If the Safety Valve should leak continually with a pressure build-up or should it cause an interruption on a sterilizing cycle prematurely (*below 124° Centigrade on the temperature gauge*), it must be replaced.


### **THE FLUE:**

The Flue serves as a protective shield for the safety valve, exhaust valve, and electrical components as well as a mounting base for the control panel. The Flue cover may be removed to allow more room for servicing the control components.

**THE EXHAUST SOLENOID VALVE:** The exhaust solenoid is normally closed. It opens at the start of heating cycle to allow cold air to escape the chamber. **It closes when chamber temperature reaches 209°F (98°C).** It also opens during a FAST VENT cycle.



## SECTION 6 ILLUSTRATED PARTS



Vent Piping View


ITEM	PART NO.	DESCRIPTION	QTY.
1	10-0938	Exhaust Solenoid Valve, 220/240 Volts	1
2	95-3730	90° Street Elbow, Brass, ½" IPS, (reworked)	1
3	10-7942	Safety Valve, 17 lbs	1
4		Heat Sink	1
5	10-2863	90° Street Elbow, Brass, 1/2" NPT	1
6	10-3327	1/2" IPS Plug, Square Head, Brass	1
7	95-6318	Compression Tube Fitting, 1/8" x 1/4", Brass (Drilled out with #29 – 0.136" Bit)	1
8		Temperature Probe (Not Shown)	1
9	15-7212	Hose, S/S Teflon, 9" Length	1
10	10-3945	Hose Clamp	2
11	95-6319	Compression Fitting, 1/8"-27 NPT, Brass (Drilled out orifice with 0.0625" bit)	1
12	10-7988	1/4 NPT- 1/8 FPT Bushing, Brass	1
13	10-3741	3/8 – 1/4 Hex Bushing	1
14	10-1054	3/8 IPS Elbow, Brass	2
15	10-1056	1/2 - 1/2 - 3/8, Tee, Brass	1
16	08-5007	1/2 - 1/8 NPT Reducer Bushing, Brass	1
17	10-3352	1/2" NPT Tee, Brass	1
*18	08-4999	1/4" NPT Plug, Brass	1
19	10-3644	3/8" IPS Plug, Square Head, Brass	1
20	10-1055	3/8 – 3/8 – 3/8 FPT Tee, Brass	1
21	95-3216	5/8" OD Copper Tubing, Soft, 9" Length	1
22	10-7987	3/8" IPS Close Nipple	1
23	10-1057	3/8" IPS Straight Union, Brass	1

\* Reserved for secondary temperature probe used only for testing calibration of autoclave.

## SECTION 6 ILLUSTRATED PARTS






*Electrical Components, Open Top View, Rear*



*Electrical Components, Open Top View, Front*

ITEM	PART NO.	DESCRIPTION	QTY.
1	09-6484	Contactora, 240V, 75 Amp	1
2		Solid State Relay, 50 Amp	1
3		Contactora, 2 Pole	3
4		Axial Fan, 230V	1
5		Printer	1
6		Controller Board	1
7		Low Water LED	1
8	10-5990	Low Water Cutoff	1
9		Power Switch	1
10		Power Supply	1
11		Fuse, Slo-Blo	1
12		Fuse Holder	1

## SECTION 6 ILLUSTRATED PARTS



*Sterilizer Assembly*

ITEM	ALL MODELS		DESCRIPTION
	QTY.	PART NO.	
1	1	95-2558	Flue Cover Assy.
2	1	95-2652	Flue Outer Case Wrap
3	1	95-2650	Upper Case, Front
4	1	95-2616	Lower Case, Front
5	1	10-6363	Insulation, Body
6	1	10-6365	Insulation, Back
7	1	10-6364	Insulation, Bottom
8	1	95-0465	Bottom Cover for Elements
9	1	95-2628	<i>Cylinder, 208V - 240V (Shown with Door Assy.)</i>





**SECTION 5 FIELD SERVICE INSTRUCTIONS & ASSEMBLY**



*Pan Supports and Baffle*

ITEM	ALL MODELS		DESCRIPTION
	QTY.	PART NO.	
1	1	95-3196	Outside Case, Left Side
2	1	95-3195	Outside Case, Right Side
3	1	95-3194	Outside Case, Back
4	1	10-0226	Handle Bumper
5	1	95-3484	Terminal Box Cover
6	2	95-2545	Pan Rack, 1 Left & 1 Right
7	1	95-2637	Condensate Baffle, Upper
8	1	95-3207	Perforated Water ( <i>Splash</i> ) Baffle
9	4	95-3284	Wear Strip

**SECTION 6 ILLUSTRATED PARTS**



## **SECTION 6 ILLUSTRATED PARTS**

### **DOOR HANDLE ASSEMBLY**

<b>ITEM</b>	<b>PART NO.</b>	<b>DESCRIPTION</b>
1	10-2318	10-32 Acorn Nut
2	10-2514	#10 Shakeproof Lockwasher
3	95-0571	10-32 Machine Screw 1 3/8" Lg.
4	95-0120	Bearing Spacer
5	95-0136	Door Lock Casting
6	10-2517	3/8" Shakeproof Lockwasher
7	10-0050	Door Lock Knob
8	95-0134	Door Handle Casting
9	10-2359	1/4"-20 Acorn Nut
10	95-0658	Door Handle Bearing Stud
11	95-0659	Door Handle Bearing Plate
12	10-2513	1/4" Shakeproof Lockwasher
13	10-1731	1/4"-20 Machine Screw 5/8" Lg.
14	95-0190	Door Lock Casting Assy. ( <i>Items 1 through 6</i> )
15	95-0145	Door Lock Knob Assy. ( <i>Items 1 through 7</i> )
16	95-0144	Complete Door Handle Assy. ( <i>Items 1 through 13</i> )
17	95-0198	Handle Bushing ( <i>Not Shown</i> )

## SECTION 7 TROUBLE-SHOOTING

### TROUBLE-SHOOTING GUIDE

TROUBLE	POSSIBLE CAUSE	CORRECTION
Sterilizer fails to operate at all ( <i>no pressure build up</i> ).	<ol style="list-style-type: none"> <li>1. Not installed correctly.</li> <li>2. Blown fuse.</li> <li>3. Contactor burned out.</li> <li>4. Wiring is defective.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check wire diagram for correct hook up.</li> <li>2. Replace fuse. If it blows, check that source of electric supply is 60 amps.</li> <li>3. Replace.</li> <li>4. Check all wiring. Repair or replace.</li> </ol>
Sterilizer operates, but fails to build up 15.5 PSI pressure.	<ol style="list-style-type: none"> <li>1. Current not heating all of the elements.</li> <li>2. Exhaust valve fails to hold pressure at 15.5 PSI.</li> <li>3. Steam leaks around door.</li> <li>4. Safety valve blows-off prematurely.</li> </ol>	<ol style="list-style-type: none"> <li>1. Remove lower front panel and see if the heating elements are working..</li> <li>2. Replace exhaust valve.</li> <li>3. Check for worn gasket or make door adjustment.</li> <li>4. Replace safety valve.</li> </ol>
Unit releases pressure before cycle has terminated on timer.	<ol style="list-style-type: none"> <li>1. Low water cut-off has functioned prematurely.</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace low water cut-off.</li> </ol>
Left and/or right side(s) heating element(s) remain on during TIMING/STERILIZATION cycle.	<ol style="list-style-type: none"> <li>1. Contactors of the temperature control switch remains closed.</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace switches.</li> </ol>



## **SECTION 8 APPENDIX**

### **A. Error Codes**

When an error occurs such as a Low Water condition, an error number will be displayed on the Digital LCD Display. Following is a list of the error numbers and their descriptions;

- **Err 01**
  - Factory calibration corrupted. Must perform factory calibration.
  - This would typically be done initially at the factory only
  
- **Err 02**
  - User setup data corrupted, user setup will be reset to defaults on any key press
  
- **Err 03**
  - Thermistors probe input open. This indicates one of three conditions;
    - The user did not set a temperature before starting a cycle
    - One or more probe wires are not connected
    - The probe has failed and should be replaced
  - Once the problem has been resolved any key press will clear the error and stop the buzzer
  
- **Err 04**
  - Thermistor probe over max temperature limit of 350°F (177°C)
  
- **Err 05**
  - PCB ambient sensor senses temperature above limit of 115°F (46.1°C)
  
- **Err 06**
  - Low water warning. Water must be added then hit the RESET button.

## **SECTION 9 WARRANTY INFORMATION**

### **STERILIZER (AUTOCLAVE) WARRANTY**

#### **MODELS: STM-E, STM-ED, STM-EL, STM-EX\*, STM-EDX\*, and STM-ELX\***

We warrant to the original purchaser that the sterilizers manufactured by Market Forge Industries, Inc. will be free from defects in material and factory workmanship if properly installed and operated under normal conditions. Within one year from date of original installation, or within 15 months from date of shipment from factory, whichever is sooner, we will repair or replace that part of any such machine that becomes defective at no cost to the customer.

This warranty is effective for One (1) Year Parts and 90 Days Labor, Travel and Mileage.

This warranty does not apply to damage resulting from use of hard or corrosive water, from failure to drain and dry cylinder daily or from inadequate cleaning procedures. Nor does it cover any part or assembly, which has been subjected to accident, alteration, or is from a machine where the serial number has been removed or altered. Normal service adjustments are not covered by this warranty.

Any defect during the warranty period shall be brought to the attention of a factory authorized service agency or the dealer from whom the equipment was purchased. He will be authorized to furnish or arrange for repairs or replacements within the terms of the warranty.

PLEASE NOTE: This warranty only applies to the USA and Canada. Elsewhere, warranty covers parts only for one year as described above.

\* *Export Model.*



Sterilizers.com - Your Sterilizer Expert  
Distribution - Maintenance - Guaranteed  
Alfa Medical - 10 Bond St Great Neck NY 11021  
1-800-762-1586 - info@Sterilizers.com



**SECTION 9 WARRANTY INFORMATION**

**STERILIZER (AUTOCLAVE) WARRANTY**

- DOMESTIC MODELS:**  **STM-E** (95-2678)  
 **STM-EL** (95-3441)  
 **STM-ED** (95-6300)

- EXPORT MODELS:**  **STM-EX** (95-2902)  
 **STM-ELX** (95-2903)  
 **STM-EDX** (95-6301)

We warrant to the original purchaser that the sterilizers manufactured by Market Forge Industries, Inc. will be free from defects in material and factory workmanship if properly installed and operated under normal conditions. Within one year from date of original installation, or within 15 months from date of shipment from factory, whichever is sooner, we will repair or replace that part of any such machine that becomes defective at no cost to the customer.

This warranty is effective for One (1) Year Parts and 90 Days Labor, Travel and Mileage.

This warranty does not apply to damage resulting from use of hard or corrosive water, from failure to drain and dry cylinder daily or from inadequate cleaning procedures. Nor does it cover any part or assembly, which has been subjected to accident, alteration, or is from a machine where the serial number has been removed or altered. Normal service adjustments are not covered by this warranty.

Any defect during the warranty period shall be brought to the attention of a factory authorized service agency or the dealer from whom the equipment was purchased. They will arrange for repairs or replacement parts within the terms of the warranty.

***PLEASE NOTE:*** *This warranty only applies to the USA and Canada. Elsewhere, warranty covers parts only for one year as described above.*

*Tear Off Bottom Section and Return to Market Forge Ind., Inc. Keep Top Section. All Fields Required.*

**Sterilizer Warranty Registration**

Model No. \_\_\_\_\_ Serial No. \_\_\_\_\_

Name \_\_\_\_\_

Company/Organization \_\_\_\_\_

Telephone No. \_\_\_\_\_ Date Installed \_\_\_\_ / \_\_\_\_ / \_\_\_\_

Street Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Email Address \_\_\_\_\_

**Mail to:** Alfa Medical - 10 Bond St Great Neck NY 11021  
1-800-762-1586 - info@Sterilizers.com

## SECTION 9 WARRANTY INFORMATION

**Disclaimer:**

*Thank you for Submitting your Sterilizer Warranty Information. This info is just for our records. Your information will not be sold. We may share your information with your local Market Forge Sales Representatives for information to respond to inquiries, inform you about products, services or events that are relevant to your sterilizer purchase. If you have any questions please don't hesitate to contact us at*

**Sterilizers.com - Your Sterilizer Expert - Distribution - Maintenance - Guaranteed**  
**Alfa Medical - 10 Bond St Great Neck NY 11021**  
**1-800-762-1586 - info@Sterilizers.com**



The serial number and model number of your unit are located on the top (roof) of the sterilizer (show in the picture on the left).

**Serial Number**

**Model Number**

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

**Sterilizers.com - Your Sterilizer Expert - Distribution - Maintenance - Guaranteed**  
**Alfa Medical - 10 Bond St Great Neck NY 11021**  
**1-800-762-1586 - info@Sterilizers.com**

